

Japet.^{w+} Technical specifications.



Technical data.

8

hours

Autonomy

of the battery under 20°C
and at 20% humidity

3

years

Lifetime

for components
of the active system

2

years

Warranty

on the active system
accessories excluded

4

kilograms

Of thrust

for each of the
4 actuators

100

%

of the active system

Assembled

internally by our
Japet experts

4

Sizes of

Second skin

available for the 2 sizes
of the active system

IIa

Classe

Device

regulation of
medical devices

5

patents

Of invention

Registered IP: FR, US
Canada, Japan, South
Korea

6

years

innovation

in research and
development

90

%

of our components


Sourced in the EU

suitcase in Shanghai



Sizes.

Hip size circumference (HC)
(On the iliac crests) **cm**




Second skin

S	UCC HC	80-89 max 86-96 max	L	UCC HC	98-110 max 106-120 max		Size 1.
M	UCC HC	89-98 max 96-106 max	XL	UCC HC	110-122 max 120-134 max		Size 2.

Active system

Schemas.



Contraindications •

- Neurological root and spinal cord deficits
- Treatment with implanted neurostimulation
- Heart or circulatory disease and severe respiratory problems
- Tumour or infectious process of the vertebrae
- Osteoporosis (< -2.5 ds)
- History of arthrodesis surgery
- Herniated disc surgery (< 3 months)
- Fracture of the dorsolumbar spine (< 3 months)
- Rib fracture (< 3 months)
- Skin injury, bruising and stretching injuries of the trunk
- Pregnancy (> 4 months)
- Previous lumbar prosthesis surgery

Textile.

Active system

Origin : **European, OEKO-TEX certified**

Velcro Lifetime : **10,000 opening and closing cycles**

Second skin

Origin : **European, OEKO-TEX certified**

Velcro Lifetime : **10,000 opening and closing cycles**

Washing recommendations : **40°, tumble drying, bleaching and ironing prohibited.**

Washing cycles : **30 valid washings**

Japet.W⁺

en.www.japet.eu

12, avenue Pierre Mauroy 59120 Loos - France

Japet.W+ is indicated for people suffering from chronic and acute low back pain. The device is intended for ambulatory distraction of the lumbar spine. Japet.W+ is a class IIa medical device manufactured by Japet Medical Devices. The CE conformity assessment has been performed by BSI Group United Kingdom (CE 2797). Read the user guide carefully.
@2024 - Japet Medical Devices - All rights reserved - Technical specifications - V.02

